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Dear Madam, Dear Sir, Dear current or prospective LiMAx user,

We are pleased to present the next edition of our quarterly LiMAx literature service today.

Each quarter, we select one or more publications from international journals, which may be of interest to you regarding our unique LiMAx method.

Have you missed a preceding edition? No problem! From now on, all editions of our literature service can be easily found on our homepage using the following link [Humedics - Literature Service](#).

This time the following publication was selected:

¹³C-methacetin breath test predicts survival in patients with hepatocellular carcinoma undergoing transarterial chemoembolization
Gairing et al., 2022. Clin Transl Gastroenterol. 2022, Sep 22. doi: 10.14309/ctg.000000000000529. Online ahead of print.

Hepatocellular cancer (HCC) – a tumor ranked third in regards of mortality among all malignant tumors worldwide - is the most common primary cancer. As the global incidence of HCC is expected to rise due to the global increase of non-alcoholic fatty liver disease (NAFLD), HCC treatments receive enhanced clinical interest. Currently, the transarterial chemoembolization (TACE) is the standard-of-care treatment of intermediate-stage HCC. Unfortunately, although this procedure was designed to minimize the harm for healthy liver tissue, TACE can reduce residual liver function in patients with HCC. As poor residual liver functions can prevent further sequential therapies, it is critical to assess function liver reserve before TACE as accurately as possible. As the prognostic power of current blood-based function parameters and derived scores for detecting TACE-based alterations in liver function is still questionable, the aim of this proof-of-concept study was to investigate the prognostic potential of the ^{13}C -methacetin breath test (^{13}C -MBT) in patients with HCC undergoing a TACE.

In total, 25 patients with a histologically and radiologically confirmed HCC who had undergone a TACE between 2017 and 2021 at the Medical University Center of the Johannes Gutenberg-University were investigated. In addition to albumin-bilirubin (ALBI), Child-Pugh and Model for End Stage Liver Disease (MELD) scores, the ^{13}C -methacetin breath test values were estimated before (pre-TACE) and after (post-TACE) the transarterial chemoembolization. To estimate the prognostic value of these parameters on the 90-day mortality, receiver operation curves (ROCs) and areas under the ROC (AUROCs) were determined.

The authors found the following:

- Pre-TACE ^{13}C -MBT values negatively correlated with ALBI and Child-Pugh scores. In addition, a strong trend - almost reached the level of significance - was found for the correlation with pre-TACE MELD scores.
- While a significant longitudinal effect (pre-TACE vs. post-TACE) was found for ALBI and MELD score, there was no time-related difference for ^{13}C -MBT values and Child-Pugh scores.
- Pre-TACE values as well as the longitudinal changes of ^{13}C -MBT, ALBI and Child-Pugh scores were significantly associated with the overall survival. For the MELD score, the longitudinal difference, but not the pre-TACE values, were significantly associated with the overall survival.
- The longitudinal change of the ^{13}C -MBT values exhibited higher 90-day mortality AUROCs, compared to the ALBI, Child-Pugh and MELD scores.
- Adverse events were not observed for the ^{13}C -MBT test.

The authors conclude:

- The study confirmed that both the baseline and longitudinal ^{13}C -MBT values predict the outcome of HCC patients undergoing TACE.
- The relative change in ^{13}C -MBT values after the first TACE may assist in identifying patients who will not benefit from further TACE treatment. It can help to identify patients in whom TACE should be discontinued because they will not derive clinical benefits from further TACE treatment.
- The lack of significant longitudinal changes in the ^{13}C -MBT values as well as correlations between ^{13}C -MBT values and ALBI scores are in line with previous investigations.
- Due to the small sample size, test results must be interpreted with caution. Especially, as the ^{13}C -MBT values were collected after the first TACE session only, the value of the ^{13}C -MBT in the context of longitudinal liver function monitoring during TACE might be underestimated in the current study.
- Due to the high clinical relevance of the survival prediction for the TACE treatment, larger clinical trials focusing on the ^{13}C -MBT value in this context should be executed.

Opinion of Humedics:

- Like previous pilot studies, this study investigated the potential of the ^{13}C -MBT (LiMax test) in predicting short-term survivability in HCC patients undergoing TACE.
- The study results showed once again that the LiMax test is a safe

procedure for estimating the liver function. They further showed the correlation between LiMAX test related values and common blood analysis parameters.

- The univariate Cox regressions as well as the AUROCs showed a better survival prediction for the ¹³C-MBT value compared to common parameters such as the MELD score. This finding underpins our previous statement (i.e. Newsletter Q1 2022) that if larger scale studies can confirm these results, the LiMAX test might be a good candidate for becoming an essential part of a future gold standard for prediction of TACE outcomes.

Please find [here](#) the abstract to this publication.

A full text copy is available on request by our [LiMAX Customer Care](#). Please get in touch with us!

If you no longer wish to receive this literature service (to: veselina.zhivkova@humedics.de), please unsubscribe [here](#).

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Humedics GmbH