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Dear Madam, Dear Sir, Dear current or prospective LiMAX user,

We are pleased to present the next edition of our LiMAX literature service today.

We regularly select one or more publications from international journals, which may be of interest to you regarding our unique LiMAX method.

Have you missed a preceding edition? No problem! From now on, all editions of our literature service can be easily found on our homepage using the following link [Humedics - Literature Service](#).

This time the following publication was selected:

The Liver Maximum Capacity Test (LiMAX) Is Associated with Short-Term Survival in Patients with Early Stage HCC Undergoing Transarterial Treatment

Fischer et al., 2022. Cancers (Basel), 14 (21):5323, DOI: 10.3390/cancers14215323

Hepatocellular cancer (HCC) is the sixth most frequent cancer and the third leading cause of cancer-related deaths worldwide. Although treatment strategies such as transarterial chemoembolization (TACE) and transarterial radioembolization (TARE) can help control tumor growth and bridge the time to liver transplantation, individual treatment response largely depends on individual factors like liver function. Subsequently, as outcome prediction is a large medical need in the context of personalized medicine, several algorithms considering parameters such as serum

bilirubin levels, Child-Turcotte-Pugh (CTP) stages, model for end-stage liver disease (MELD) and albumin-bilirubin (ALBI) have been developed in the past. Nevertheless, as these algorithms depend on laboratory parameters and the subjective estimation of clinical symptoms, the aim of this study was to determine the correlation of a direct liver function measurement (LiMAX test) with short-term survival in patients with early and intermediate stage HCC after TACE and TARE.

In total, 69 patients who underwent a TACE or TARE at the University Medical Center Leipzig between November 2017 and April 2020 were investigated. Laboratory data and cancer characteristics, CTP-, MELD and ALBI scores as well as the LiMAX test values and well-established serological analyses were determined before and 4 weeks after the TACE or TARE procedures. In addition, the occurrence of adverse events at 4 weeks after TACE/TARE and survival until 34 months after treatment were assessed. Lastly, beside common inferential statistics, survival analyses were performed with Kaplan-Meier estimators and Cox regressions analysis for the 30-, 60-, 90- and 180-day survival.

The authors found the following:

- Pre-treatment LiMAX values significantly correlated with bilirubin, albumin, aspartate aminotransferase (AST), international normalized ratio (INR) as well as ALBI and MELD scores. After treatment, these correlations remained for bilirubin, albumin as well as ALBI and MELD scores. In addition, a correlation of post-treatment LiMAX values and CTP occurred.
- LiMAX values, blood parameters as well as the MELD, ALBI and CTP scores did not significantly change between the day before and 4 weeks after transarterial treatment.
- Significantly different LiMAX values were observed in subgroups stratified by ALBI and CTP grades.
- Patients with LiMAX values ≤ 150 $\mu\text{kg/h}$ showed significantly increased levels of liver enzymes, ALBI, MELD and CTP scores as well as decreased levels of albumin and platelet counts compared with patients with higher LiMAX (>150 $\mu\text{kg/h}$) values before and after treatment. They further showed significantly lower 30- and 60-day survival. Despite this, the differences between these groups for the 90- and 180-day survival as well the overall survival were indistinct.
- In patients with BCLC stage A, but not for stage B and C, with LiMAX values ≤ 150 $\mu\text{kg/h}$ showed significantly lower survival rates compared to patients with higher LiMAX (>150 $\mu\text{kg/h}$) values.
- In univariate Cox regression analysis, the overall survival was associated with pre-treatment bilirubin and albumin levels, CTP Child B as well as MELD, ALBI score grade 2 and 3. Although MELD and ALBI grade 3 remained predictors for the overall survival in multivariate Cox regression analysis, after Bonferroni correction for multiple testing, none of the parameters was associated with overall survival.
- No significant association between the occurrence of common adverse events and pre-treatment LiMAX values was found.
- The incidence of complication was associated with pre-treatment CTP and ALBI scores as well as AST levels in univariate but not in multivariate regression analysis.

The authors conclude:

- As LiMAX values ≤ 150 $\mu\text{kg/h}$ are strongly associated with decreased survival rates over 30 and 60 days as well as over 180 days in patients with BCLC stage A, LiMAX measurements might be a useful tool for identifying patients with short-term survival.
- Due to the relationship of lower LiMAX values and BCLC stages on long-term survival, it needs to be investigated whether those patients might have drawn more benefit from different treatment approaches.
- As severe complications of TACE such as liver failure, biliary or hepatic artery injury and infection as well as mortality are significantly less than common side effects such as the post-embolization syndrome, this study was not able to reliably detect associations between the LiMAX test and such rare events.
- Although several studies showed that the LiMAX test was appropriate to quantify the liver function capacity in different clinical settings, and although these studies also exhibited the correlation of pre-treatment LiMAX values with surrogate markers indicating post treatment liver function deterioration, this study found no differences between LiMAX values before and 4 weeks after TACE or TARE. Hence, it seems that the LiMAX test might only be successful to detect short-term changes in liver function and will not work for time periods where the patients fully recovered after treatment.

Opinion of Humedics:

- This study showed the potential of the LiMAx test in the field of HCC treatment. Once again, it has been shown that LiMAx values correlate with common outcome predictors such as bilirubin, albumin, aspartate aminotransferase (AST), international normalized ratio (INR) as well as ALBI and MELD scores.
- The significant lower 30 to 60-day survival for patients with lower LiMAx values shows the potential of the LiMAx test for survival prediction.
- Due to the small number of severe complications, the sample size of this study was not able to detect the prediction potential of the LiMAx test on the occurrence of severe complications such as liver failure, biliary or hepatic artery injury or infections. We would be pleased if larger scale studies would investigate the potential of the LiMAx test on the prediction of severe complications in patients with early stage HCC undergoing a transarterial treatment.

For a summary review of the scientific data on the LiMAx test, we recommend the presentation by Dr. Kohnen at this year's **Paris Hepatology Conference (PHC)**. The presentation can be found on our homepage via [Humedics - Events](#).

The Humedics team wishes you a happy Easter.

Please find [here](#) the abstract to this publication.

A full text copy is available on request by our [LiMAx Customer Care](#). Please get in touch with us!

If you no longer wish to receive this literature service (to: benny.maulana@humedics.de), please unsubscribe [here](#).

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Humedics GmbH